## Name: General Use Case

Actor(s): User, HeartWave (Device)

Precondition(s): None

Success Guarantee: User is able to start a session where the device rgiecords the user's HRV and possibly other metrics. User is then able to end the session whenever they want and get a detailed summary of the info recorded by the device.

## Main success scenario:

- 1. User attaches the sensor to their body
- 2. User presses the on button on the device
- 3. Device turns on and displays a menu on the screen
- 4. User selects the Challenge option
- 5. Device screen displays challenge levels and breath pace settings
- 6. User selects a challenge level
- 7. Device screen displays challenge info and begins capturing user's body data
- 8. User manually ends the session
- 9. Device screen displays summary of session info

## Extensions:

6a. User can select a breath pace level (between 1-30 seconds default 10s)

7a. If the breath pacer is attaches, device also displays breath data

## Misc Use Cases that I'm too lazy to format:

- 1. Other menu actions: (View session log/history, starting a new session, go to settings)
- 2. Reset/wipe data
- 3. Low battery